



**Training Manual
on
Water Quality Testing &
Report Management (WQTRM) Software**



**MANUAL DEVELOPED BY:
ENGR. MD. JAHIDUL HASAN, CONSULTANT**

Office of the Executive Engineer
Training Department
HRD Building
36-37, Bir Uttam AK Khandkar Road, Mohakhali, Dhaka-1212

Manual Compilation

Engr. Md. Jahidul Hasan

Printing Support :

Emergency Multi-Sector Project for Rohingya Crisis Response (EMCRP)

Coordinator:

Training Department of Public Health Engineering (DPHE), Mohakhali, Dhaka-1212

Publication Date :

February 2024

Design & Printing :

Emergency Multi-Sector Project for Rohingya Crisis Response (EMCRP)

Copywrite:

Department of public health engineering

TABLE OF CONTENTS

SESSION	CONTENT	PAGE NO
First session:	<ul style="list-style-type: none"> ● Introduction and purpose of training ● Ceation of learning environment. 	
Second session:	<ul style="list-style-type: none"> ● Fundamental concept of Computer Hardware, Software ● Computer Function ● Computer System 	
Third session:	<ul style="list-style-type: none"> ● What is Water Quality Testing and Report Management (WQTRM) Software? ● Current limitations of WQTRM Software, ● Outcome or benefit of the new System, ● Main Features of the Proposed System. 	
Fourth session:	<ul style="list-style-type: none"> ● Operational flow chart of (WQTRM) System ● What are Software Modules & Main Modules of WQTRM ● Technical features and system Architect of WQTRM Software. 	
Fifth session:	<ul style="list-style-type: none"> ● Master data & Master Data entry ● User Management ● Role Management ● Address Management. Responsibilities of System Administrator 	
Sixth session:	<ul style="list-style-type: none"> ● What is Transecting Data Entry ● Customer Registration Entry ● Application Management, Sample Collection ● Bill Management ● Responsibilities of System SAE and EE 	
Seventh sessions:	<ul style="list-style-type: none"> ● Lab Management System (LIS), ● Responsibilities of LAB Technician, Chemist LAB In-Charge 	
Eight sessions:	<ul style="list-style-type: none"> ● Payment Gateway Integration ● Mobile Base Application ● Notification Management System ● Dynamic Report ● Dynamic Dashboard. 	

Foreword

Immediate provision of emergency services like water supply, sanitation, and health services in case of any disaster or emergency is a very challenging task. The Rohingya crisis in Cox's Bazar is no exception. Providing water supply, providing sanitation facilities, proper waste management, and promoting hygiene in the densely populated Rohingya camps is a huge challenge. Under the pressure of this huge population, Ukhia and Teknaf upazilas are facing a serious crisis in the quality of life along with the natural, socio-economic environment and overall development of Cox's Bazar district is being hindered. Therefore, considering the overall situation of Cox's Bazar district, with the financial support of the World Bank, the Directorate of Public Health Engineering is implementing the "Multi-Sector Rohingya Crisis Response" project in the entire Cox's Bazar district, including Ukhia and Teknaf Upazilas, which provide Rohingya camps and shelters. Under this project, the Department of Public Health Engineering is conducting multi-faceted activities including climate-resilient safe water supply, sanitation, sewage management and solid waste management, and promotion of .social and gender-based integrated hygiene

Nine training manuals have been formulated to increase the knowledge and skills of the officers, employees, and related parties of the Directorate of Public Health Engineering, Project Office, Consultant, Implementation and Management and Maintenance Institutions associated with the implementation of the project activities to achieve the desired goals through the proper management of the above activities. Skill Development Manual on " Water Quality .Testing & Report Management (WQTRM) System" is part of it

Through this manual, the trainees will learn the correct information about the various components and mechanisms of using solar energy. The training course is specifically aimed at .making learning interesting, interactive, and practical

Sincere thanks to the Hon'ble Chief Engineer, Department of Public Health Engineering, and Project Director, EMCRP for providing comprehensive and comprehensive guidance to the Training Division, DPHE in preparation of this manual. Md. Jahidul Hasan, short-term consultant, EMCRP and Md. Abdullah-Hil-Kafi, Assistant Engineer, Training Division, DPHE, has directly or indirectly contributed to the preparation of this training manual. I hope that this manual will not only be a single topic related to the EMCRP project but will contribute as a .complete guideline for DPHE's solar energy utilization work in the future



Ruksana Parvin
Executive Engineer
Training Division
Department of Public Health Engineering (DPHE)

Sequential

Since August 2017, a large part of the Rohingya population forcibly displaced from Myanmar has entered the coastal district of Cox's Bazar in Bangladesh, creating one of the world's largest displacement crises. About 1.1 million displaced Rohingya have settled in makeshift camps in Ukhia and Teknaf upazilas—almost three times the local population. The area has become very densely populated due to the shelter of the displaced Rohingya population. Besides, the infrastructure there is very weak, lack of availability of basic services, and prone to a high risk of natural disasters, especially cyclones and floods

Analyzing the initial humanitarian assistance operations, the World Bank has coordinated its aid-intensive ongoing operations with long-term operations for the displaced Rohingya population. As part of this, the "Emergency Multi-Sector Project for Rohingya Crisis Response (EMCRP)" is currently being implemented by the Directorate of Public Health Engineering to provide safe water, sanitation, and hygiene facilities to the displaced Rohingya population. The main objective of this project is to provide safe water supply and sanitation services as well as promote social and gender-integrated sanitation

To increase the knowledge and skills of the officials, consultants, implementers, and management and maintenance institutions at various levels of the Department of Public Health Engineering, related to the implementation and management of proper safe water supply, sanitation services, sanitation and sewage management, etc. Nine training manuals have been prepared on relevant subjects. I hope that these manuals will make a significant contribution to improving the skills of all concerned

I express my gratitude to all the persons who contributed directly or indirectly by providing valuable feedback and suggestions for the preparation of the training manuals. I am grateful to the World Bank for providing grant-in-aid support to the Multi-Sector Project on Rohingya Crisis Mokobila on an urgent basis. Congratulations and thanks to all the colleagues the Executive Engineer, the Training Department, the short-term Consultant, EMCRP, and the Project Director of Public Health Engineering for their cooperation in finalizing and preparing these manuals



Md. Sarwar Hossain

Chief Engineer

Department of Public Health Engineering (DPHE)

Acknowledgment

Since August 2017, Rohingyas have crossed the border into Bangladesh's Cox's Bazar district due to extreme violence in Myanmar's Rakhine State. The Bangladesh government provided them temporary shelter in Ukha and Teknaf upazilas of Cox's Bazar district on humanitarian grounds. As a result of the location of this large number of displaced Rohingya people, the Rohingya camp area of Cox's Bazar has become the most densely populated area in the world, and meeting the basic needs of the population has become a big challenge. As a result, the development of Cox's Bazar district is hindered and the quality of life of the local people is threatened. To deal with this crisis, the "multi-sector project to deal with the Rohingya crisis on an emergency basis" is being implemented by the Department of Public Health Engineering on behalf of the Government of Bangladesh to assist with safe water supply, sanitation, and health development. The Directorate of Public Health Engineering EMCRP and the Training Department are working together to provide the necessary support to increase the skills of the concerned officers/employees under the department and various workers under the scope of the project

As a part of this work, nine training manuals have been formulated on various relevant topics to increase the knowledge and skills of the officials, employees, office assistants, operators, and supervisors of the Department of Public Health Engineering, Project Office, Consultants, Implementation and Management and Maintenance Institutions associated with the implementation of EMCRP project activities. Skill Development Manual on " Water Quality Testing & Report Management (WQTRM) System " is part of it. This training manual includes a detailed discussion on the relevant topics, session planning, and various required materials (learning and reference materials/reading materials) considering the participants. I strongly believe that through the manual, the activities of Water Quality utilization under the project will be more standardized and the project activities will be accelerated

I would like to express my gratitude to all those who contributed their valuable time, feedback, and suggestions including Ruksana Parvin, Executive Engineer, Training Division, DPHE, and Md. Jahidul Hasan, Short Term Consultant, EMCRP. I express my gratitude to Mr. Md. Sarwar Hossain, Honorable Chief Engineer, Department of Public Health Engineering, for his utmost cooperation in the preparation of the manual by providing overall guidance. I am grateful to the World Bank authorities for providing overall support in project management through principal and additional financing of the project. Also, special thanks to EMCRP Project Social Development Officer Mr. Md. Muktedir Harun, Deputy Project Director and Project Director's office consultants, training consultant Mr. Md. Shahidur Rahman, and related officials for their overall cooperation in finalizing the manual. I hope that through this training manual, those involved in the said work will be properly informed and acquire the necessary skills



Mohammad Abdul Qayyum

Project Director (EMCRP)

Department of Public Health Engineering (DPHE)

TRAINING SCHEDULE

1st Day (8 Hours)

TIME	SEASON	TOPIC
09.00 - 09.30	Registration	
09.30 - 10.00	First session 30min	<ul style="list-style-type: none"> ● Introduction and purpose of training ● Creation of learning environment.
10.10 - 10.30	Tea break 30min	
10.30 - 11.30	Second session 1hr 30 min	<ul style="list-style-type: none"> ● Computer, Hardware, Software ● Computer Function, Computer System
11.30 - 12.30	Third session 1hr	<ul style="list-style-type: none"> ● What is Water Quality Testing and Report Management (WQTRM) Software ● Current limitations of WQTRM Software: ● Outcome or benefit of the new system ● Main Features of the Proposed System
12.30 - 1.30	Fourth session 1hr	<ul style="list-style-type: none"> ● Operational flow chart of (WQTRM) System ● What are Software Modules & Main Modules of WQTRM ● Technical features and system Architect of WQTRM Software
01.30 - 02.30	Food and prayer	
02.30 - 04.30	Fifth session 1.30hr	<ul style="list-style-type: none"> ● Master data & Master Data entry ● User Management ● Role Management ● Address Management ● Responsibilities of System Administrator
04.30 - 05.00	Tea break	

First Day



First session

Describe the purpose of the training
Pre-training assessment

First session

1.1 Introductory episode:

- ◆ Describe the purpose of the training
- ◆ Pre-training assessment

1.2 Introductory Session Objectives: At the end of this session the participants will-

- ◆ Participants will get to know each other
- ◆ Get an idea about the objectives of the training
- ◆ A learning environment will be created
- ◆ Coordination between expectations and topics can be achieved

1.3 Duration: 30 Min

1.4 Training Materials: Whiteboard, whiteboard markers, multimedia screen, duster and manual.


1.5 Topics to be discussed:

1.5.1 Overall Objectives of Training: After this training, the participants will.

- ◆ Get an idea about Computers, Computer Hardware, Software, and Computer Systems.
- ◆ Get an idea about Water Quality Testing and Report Management (WQTRM) Software.
- ◆ Get an idea about the Present System and Proposed System.
- ◆ Get an idea about the Operational flow chart of (WQTRM) System.
- ◆ Can explain about Main Modules, Technical features, and system Architect of WQTRM Software.
- ◆ Responsibilities of Sub Assistant Engineer (SAE) and Executive (EE)
- ◆ Get an idea about the Master Data Entry System, Responsibilities of System Administrator
- ◆ Get an idea about the System (LIMS).
- ◆ Responsibilities of LAB Technician, Chemist, LAB In-Charge
- ◆ Get an idea about Discussion on the Dashboard, Reports Management System and Challenges.

1.5.2 Specific objectives of this training: At the end of this training the participants will –

- ◆ Can explain What is Computer, Computer Hardware, Software, Computer System?
- ◆ Can explain about Water Quality Testing and Report Management (WQTRM) Software.
- ◆ Can explain about the Present System and new System.
- ◆ Can tell about the Operational flow chart of (WQTRM) System.
- ◆ Can tell about Main Modules, Technical features, and system Architect of WQTRM Software.
- ◆ Can explain Master Data Entry System.
- ◆ Can explain about LAB Management system (LIMS).
- ◆ Can explain the dashboard, Reports Management System, and Challenges.



Second session

**Computer, Hardware, Software
Block Diagram, Computer Function
and Computer System**

Second session

2.1 Topic:

- ◆ Computer, Hardware, Software
- ◆ Block Diagram, Computer Function and Computer System

2.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can tell about Computer, Hardware, Software, Block Diagram and Computer Function, and Computer System

2.3 Duration: 1hr

2.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual

2.5 Procedures for Conducting Training

SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none">● The name of the subject to be discussed must be presented in writing.● What the participants want to learn from this discussion should be known through questions.● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none">● The topic should be presented through questions.● Participants should be encouraged to participate in discussions.● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none">● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none">● To assess how much the participants know the subject matter by presenting and answering questions.

2.6.1 Computer?

The computer is an advanced electronic device that

- ◆ takes raw data as input from the user
- ◆ processes these data under the control of a set of instructions (called a program)
- ◆ gives the result (output)
- ◆ and saves output for future use

A computer requires both hardware and software **Hardware:** All physical parts (that we can touch) of the computer are known as Hardware (Keyboard, Monitor, Mouse).

- ◆ Input Device: Mouse, Key Board
- ◆ Output Device: Monitor, Printer
- ◆ Storage Device: Hard Disk, Pen Drive

Software: Programs (program is a series of instructions), Software is a set of instructions, data, or programs used to operate computers and execute specific tasks. It is the opposite of hardware, which describes the physical aspects of a computer.

There are three basic types of software: Application programs, device drivers, and operating systems. Each type of software performs a completely different job, but all three work closely together to perform useful work.

Database: A database is a collection of data, typically stored in a computer system. A database is usually controlled by a database management system (DBMS).

2.6.2 Block Diagram of Computer Function-

A computer has four functions:

Input (Data): Input is the raw information entered into a computer from the input devices (Keyboard, Mouse). It is the collection of letters, Numbers, images, etc.

1.Process: Process is the operation of data as per given instruction in the CPU. It is an internal process of the computer system.

2.Output: Output is the processed data given by the output device (Printer, monitor) after data processing. Output is also called Result.


3.Storage: We can save these results in the storage devices (Hard disk, pen drive) for future use.

2.6.3 Computer/Computerized System:

All of the components of a computer system can be summarized with simple equations.

COMPUTER SYSTEM = USER + HARDWARE + SOFTWARE

- ◆ Hardware = Internal Devices + Peripheral Devices
- ◆ Software = Program
- ◆ USER = Person, who operates a computer



Picture: Block Diagram of Computer Function



Picture: Computer System

Third session

**What is Water Quality Testing and Report Management (WQTRM) Software
Outcome or benefit of the new software
Main Features of the new software**

Third session

3.1 Topic:

- ◆ What is Water Quality Testing and Report Management (WQTRM) Software
- ◆ Outcome or benefit of the new software
- ◆ Main Features of the new software

3.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can explain about Water Quality Testing and Report Management (WQTRM) Software
- ◆ Can explain about benefit of the new system
- ◆ Can Explain about new (WQTRM) Software
- ◆ Main Features of the new software

3.3 Duration: 1hr

3.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

3.5 Procedures for Conducting Training

SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none">● The name of the subject to be discussed must be presented in writing.● What the participants want to learn from this discussion should be known through questions.● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none">● The topic should be presented through questions.● Participants should be encouraged to participate in discussions.● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none">● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none">● To assess how much the participants know the subject matter by presenting and answering questions.

3.6.1 What is WQTRM Software?

WQTRM Software is one of the application software which is Developed for water quality testing all over Bangladesh and this project is funded by the World Bank.

This system will store in the central database management system application data, financial data, and LAB data for future analysis.

The main goal of **(WQTRM) Software** is to transform a conventional manual water quality testing report into a modern integrated computerized system and provide a good quality LAB report to our valuable customers.



3.6.2 Current limitations of old WQTRM Software:

- The water sample and water point are not linked
- No tracking of the sample status
- NO process of authentication/validation of water sample tracking
- Desktop-based solution
- Not integrated
- No dynamic reports
- No Dashboard

3.6.3 Outcome or benefit of the new software:

- To store application data, financial data, and Lab data in a central server for future research and real-time reporting.
- To digitalize the application process of water quality tests.
- To Computerized payment management
- To digitize sample collection, and diagnosis Reports
- To collect water point location through mobile GPS and pictures
- To provide SMS and email notifications
- Report sending from any Lab to its higher-level office
- Automated Real-time Reports management
- Real-Time Dashboard
- Automated Notifications

3.6.4 Main Features of the Proposed System

- User friendly
- Wave base solution
- Power full MIS Dashboard
- Different types of reports.
- Integrated Solution
- Central Database
- Client-server, three-tier, and cloud-based use
- High-quality user requirement-wise reports
- Better interface



Fourth session

Operational flow chart of (WQTRM) software
What is Software Modules & Main Modules of WQTRM
Technical features and system Architect of WQTRM Software

Fourth Session

4.1 Topic:

- ◆ Operational flow chart of (WQTRM) software
- ◆ What is Software Modules & Main Modules of WQTRM
- ◆ Technical features and system Architect of WQTRM Software

3.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can explain about Water Quality Testing and Report Management (WQTRM) Software
- ◆ Can explain about Main Modules of WQTRM
- ◆ Can explain about technical features and system Architect of WQTRM Software

3.3 Duration: 1hr


3.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

3.5 Procedures for Conducting Training


SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none"> ● The name of the subject to be discussed must be presented in writing. ● What the participants want to learn from this discussion should be known through questions. ● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none"> ● The topic should be presented through questions. ● Participants should be encouraged to participate in discussions. ● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none"> ● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none"> ● To assess how much the participants know the subject matter by presenting and answering questions.

4.6.1 Operational flow chart of WQTRM system

There are two types of clients- 1. Individual Clients. 2. Institutional Client.



Flow chart of Individual client



4.6.2 Software Modules?

In computer software, a module is a unit of program dedicated to a specific function.

WQTRM Software is a cost-effective, efficient software system to monitor water quality using software, hardware, and databases. It measures if water is safe to drink, fish from, swim in and if it complies with legal requirements. Water quality is measured by several factors, such as the concentration of dissolved oxygen, bacteria levels, the amount of salt (or salinity), or the amount of material suspended in the water (turbidity).

4.6.3 Technical Features of WQTRM Software

- ◆ WQTRM Software is fully scalable and adaptable to future needs, technological evolutions,
- ◆ and volumes of business
- ◆ WQTRM Software is Platform independent and should not be constrained to a single Hardware Platform Operating System or database
- ◆ WQTRM Software is Developed using open-source technology
- ◆ Recommend the detailed infrastructure with specifications required for rolling out the solution in development, testing, staging, production, high
- ◆ It is using Developing tools like JAVA, .NET, and PHP for easy and simple integration for
- ◆ third-party external systems
- ◆ Integrating with SMS and payment gateway, availability of DC/DR technology

4.6.4 Main Software Modules of WQTRM Software

1. Master Data Entry Management-

- ◆ User Management & Role Management
- ◆ Location Management
- ◆ Service Management and Rate Management

2. Transaction Entry Data

- ◆ Application Management
- ◆ Bill Management
- ◆ Lab Management
- ◆ Notification Management
- ◆ Payment Gateway Integration.
- ◆ Development of Mobile Base Application for Data Entry

3. Reports & Dashboard


- ◆ Dynamic Report System
- ◆ Dynamic Dashboard.




Main software modules

4.6.5 System Architect of WQTRM software:


System architecture is, simply, the organization of a system. This organization includes all components, how they interact with each other, the environment in which they operate, and the principles used to design the software.



System Flow Chart-01



System Flow Chart-01





Fifth session

**Master Data Entry System, User Management
Role Management, Address Management
Service Management and Rate Management
Responsibilities of System Administrator
Operation Manual of Management Admin**

Fifth Season

5.1 Topic:

Master Data Entry System

- ◆ User Management
- ◆ Role Management
- ◆ Address Management
- ◆ Service Management and Rate Management
- ◆ Responsibilities of System Administrator
- ◆ Operation Manual of Management Admin

5.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can learn about Master Data Entry System
- ◆ User Management
- ◆ Role Management
- ◆ Address Management
- ◆ Service Management and Rate Management
- ◆ Can learn responsibilities of System Administrator & Management Admin

5.3 Duration: 1hr

5.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

5.5 Procedures for Conducting Training

SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none">● The name of the subject to be discussed must be presented in writing.● What the participants want to learn from this discussion should be known through questions.● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none">● The topic should be presented through questions.● Participants should be encouraged to participate in discussions.● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none">● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none">● To assess how much the participants know the subject matter by presenting and answering questions.

5.6.1 Master Data Entry System?


Master data is the core data that is essential for running the whole computerized system. It is data about key business entities that provide context for business transactions and operations. Data designated as master data can vary from organization to organization.

Master Data Management:

1. User Management
2. Role Management
3. Address Management
4. Service Management and Rate Management

1. User Management:

User management is a system to handles activities related to individuals' access to devices, software, and services. It focuses on managing permissions for access and actions as well as monitoring usage. Functions of user management include: Providing users with authenticated access.



These are three types of Users:

System administrator- all access to the system.

admin user- of the system to create/update/delete different types of user accounts like Admin, Accounts, Laboratory, and Top Management, etc. to access and perform different activities of the system.

General User- to enter data of application, billing, Laboratory, etc.

A robust user management system is required, which will enable the authorized admin user of the system to create/update/delete different types of user accounts like Admin, Accounts, Laboratory, and Top Management, etc. to access and perform different activities of the system.

2. What are the roles of a user?

A user role is a set of privileges that is assigned to a user or user group to allow the user or user group to perform certain tasks and manage certain sessions. To be assigned to a role, each user or group of users must have a valid user ID or group ID in the user registry on the management server.

3. Home addresses Definition: Define Home addresses. means the dwelling location at which an individual resides and includes the physical address, mailing address, street address, parcel ide.

4. What is service management? : Service Management is a customer-focused approach to delivering information technology. Service Management focuses on providing value to the customer and also on the customer relationship.

1. Operation Manual of System Admin

1. To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login> and enter User Name and Password.



2. In admin dashboard, can view 'Sample Collection this Year' with sinusoidal and pie chart representations and filter for a specific date range.



3. In the "Role Manage" section, can create a role by entering a role name, selecting permission name and guard name, and submitting the form.

The 'Role Management' section shows a table with the following data:

Sl.	Name	Permissions	Action
1	system.admin	<ul style="list-style-type: none"> Sample RoleManagement SettingManagement LabList LabCreate LabDestroy LabShow LabUserCreate ReportManagement LocationManagement LabManagement LabAdmin LabUpdate LabParameter LabUser LabUserEdit 	<ul style="list-style-type: none"> + + +

A red box highlights the '+ Create' button in the top right corner.

4. In the "Employee Manage-List" section, can create and edit employee information by required data

The 'Employee Information' section shows a table with the following data:

Sl	Image	Name	Designation	Phone	Email	Signature	Chair	Action
1		Ferhan Savit	Chief Chemist	01777444501	ferhan_savit@gmail.com		<ul style="list-style-type: none"> + + 	<ul style="list-style-type: none"> + + +
2		Rashed Koirin	Chief Chemist	01777444502	rashed_koirin@gmail.com		<ul style="list-style-type: none"> + + 	<ul style="list-style-type: none"> + + +

A red box highlights the '+ Create Employee' button in the top right corner.

5. In the "Settings-DPHE" section, can create and edit circles by required data

The 'Settings-DPHE' section shows a table with the following data:

Name	ISB Code	Status	Action
Chilling Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Spill Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Whirl Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Refriger Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Refriger Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Chilling H&B Steam Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Freezing Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +
Chilled Circle	0	<input type="checkbox"/>	<ul style="list-style-type: none"> + + +

A red box highlights the '+ Create' button in the top right corner.

6. In the "Settings-Type of Water Point" section, can create and edit type of water point.



The screenshot shows a form titled "Create Dphe Circle" with a close button (X) in the top right corner. The form contains four fields: "Name" with the placeholder "Enter Type Name", "GEO Code" with the placeholder "Enter Type Others ...", "Type" with the value "Dphe Circle", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button and a red arrow pointing to it.

7. In the "Settings-Ownership Type" section, can create and edit the Ownership type.



The screenshot shows a form titled "Create Water Point Type" with a close button (X) in the top right corner. The form contains four fields: "Name" with the placeholder "Enter Type Name", "GEO Code" with the placeholder "Enter Type Others ...", "Type" with the value "Water Point Type", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button and a red arrow pointing to it.

8. In the "Settings-Area and Land Type" section, can create and edit Area and Land type.



The screenshot shows a form titled "Create Area Lend Type" with a close button (X) in the top right corner. The form contains four fields: "Name" with the placeholder "Enter Type Name", "GEO Code" with the placeholder "Enter Type Others ...", "Type" with the value "Area Lend Type", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button and a red arrow pointing to it.

9. In the "Settings-Water Source" section, can create and edit the Water Source.



The screenshot shows a form titled "Create Water Source" with a close button (X) in the top right corner. The form contains four fields: "Name" with the placeholder "Enter Type Name", "GEO Code" with the placeholder "Enter Type Others ...", "Type" with the value "Water Source", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button and a red arrow pointing to it.

10. In the "Settings-Sample Source" section, can create and edit the Sample Source.

The screenshot shows a web form titled "Create Sample Source" with a close button in the top right corner. The form contains four fields: "Name" with a placeholder "Enter Type Name", "GEO Code" with a placeholder "Enter Type Others ...", "Type" with a dropdown menu showing "Sample Source", and "Status" with a toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button. A red arrow points to the "Create" button, and a red box highlights it.

11. In the "Settings-Sample Category" section, can create and edit the Sample Category.

The screenshot shows a web form titled "Create Sample Category" with a close button in the top right corner. The form contains four fields: "Name" with a placeholder "Enter Type Name", "GEO Code" with a placeholder "Enter Type Others ...", "Type" with a dropdown menu showing "Sample Category", and "Status" with a toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button. A red arrow points to the "Create" button, and a red box highlights it.

12. In the "Settings-Bottle Type" section, can create and edit delete the Bottle Type.

The screenshot shows a web form titled "Create Bottle Type" with a close button in the top right corner. The form contains four fields: "Name" with a placeholder "Enter Type Name", "GEO Code" with a placeholder "Enter Type Others ...", "Type" with a dropdown menu showing "Bottle Type", and "Status" with a toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button. A red arrow points to the "Create" button, and a red box highlights it.

13. In the "Settings-Latrine Type" section, he can create and edit Latrine Type.

The screenshot shows a web form titled "Create Latrine Type" with a close button in the top right corner. The form contains four fields: "Name" with a placeholder "Enter Type Name", "GEO Code" with a placeholder "Enter Type Others ...", "Type" with a dropdown menu showing "Latrine Type", and "Status" with a toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button. A red arrow points to the "Create" button, and a red box highlights it.

14. In the "Settings-Latrine Platform Condition" section, can create and edit Latrine Platform Condition.

Create Latrine Platform [X]

Name * Enter Type Name

GEO Code Enter Type Others ...

Type * Latrine Platform

Status

Close Create

15. In "Settings-Install By" section, can create and edit the Install By.

Create Install By [X]

Name * Enter Type Name

GEO Code Enter Type Others ...

Type * Install By

Status

Close Create

16. In "Settings-Source Category" section, can create and edit Source Category.

Create Source Category [X]

Name * Enter Type Name

GEO Code Enter Type Others ...

Type * Source Category

Status

Close Create

17. To create a Source Category, enter required data and click "Create".

Create Source Category [X]

Name * Enter Type Name

GEO Code Enter Type Others ...

Type * Source Category

Status

Close Create

18. In "Settings-Source Sub Category" section, can create and edit Source Sub Category.




The screenshot shows a modal window titled "Create Sub Source Category" with a close button (X) in the top right corner. The form contains three fields: "Name" with a placeholder "Enter Sub Source Category Name", "Source Category" with a placeholder "Revinew Budget", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

19. In "Settings-Preservative Used" section, can create and edit the Preservative Used.



The screenshot shows a modal window titled "Create Preservative" with a close button (X) in the top right corner. The form contains four fields: "Name" with a placeholder "Enter Type Name", "GEO Code" with a placeholder "Enter Type Others ...", "Type" with a placeholder "Preservative", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

20. In "Settings-Parameter" section, can create and edit Source Sub Category.



The screenshot shows a modal window titled "Create Parameter" with a close button (X) in the top right corner. The form is divided into two columns. The left column has fields for "Name" (placeholder: "Enter Parameter Name"), "Per Unit Cost" (placeholder: "Enter Unit Price"), "Bangladesh Standard" (placeholder: "Enter Unit Price"), and "Analysis Procedure" (checkbox). The right column has fields for "Short Name" (placeholder: "Enter Parameter Short Name"), "Unit" (placeholder: "Enter Unit"), and "BDS Number" (placeholder: "Enter BDS Number"). Below these columns is a "Status" toggle switch, currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

21. In "Settings-Procedure" section, can create and edit the Procedure.



The screenshot shows a modal window titled "Create Procedure" with a close button (X) in the top right corner. The form contains three fields: "Name" with a placeholder "Enter Procedure Name", "Full Name" with a placeholder "Enter Procedure Full Name", and "Status" which is a toggle switch currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

22. In "Settings-Analysis Method" section, he can create and edit Analysis Method.

Create Analysis Method

Name *

GEO Code

Type *

Status

[Close](#) [Create](#)

23. In Lab-List , can create new lab, edit lab information, view lab parameters, view lab users, view lab details, view chemist committee, remove lab.

Lab List

ID	Name	Office/Officer	Phone	Covered Material	Address	Status	Action
1	Management Zonal Lab	Senior Chemist	0155 4333716 (No.)		301 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
2	Headquarter Zonal Lab	Senior Chemist	0227-69770		302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
3	Merut Zonal Lab	Senior Chemist	0949 297993		302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
4	Tongk Zonal Lab	Senior Chemist	02 22489007	• Stalpas	302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
5	Amambak Zonal Lab	Senior Chemist	0434 64470		302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
6	Agilok Zonal Lab	Senior Chemist	0521 63416		302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete
7	Segre Zonal Lab	Senior Chemist	031 35205		302 Labha Street, Port Saeed, RI 22480	<input checked="" type="checkbox"/>	Edit View Delete

24. To create a laboratory, enter required data and click "Save Information".

Create New Laboratory

Name Single*

Name Single*

Address

Short Name*

Phone No.*

Fax

Office Of The

[Save Information](#)

25. To create a Lab-Admin, enter required data and click "Create".

Lab Admin

SI	Name	Email	Lab	Office Of The	Employee
1	Tongk Lab Admin	ton_lab_admin@gmail.com	Tongk Zonal Lab	Senior Chemist	Employee

26. In Location-Division he can create new divisions, edit, remove divisions.



The screenshot shows a web form titled "Create a New Division" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Division Name English", "Name Bangla" with a placeholder "Enter Division Name Bangla", "GEO Code" with a placeholder "Enter District GEO Code", and "Website URL" with a placeholder "Enter Website URL". There is a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button.

27. In Location-District he can create and edit new districts.



The screenshot shows a web form titled "Create a New District" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter District Name English", "Name Bangla" with a placeholder "Enter District Name Bangla", "GEO Code" with a placeholder "Enter District GEO Code", "Division" with a placeholder "Select Division", and "Website URL" with a placeholder "Enter Website URL". There is a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button.

28. In Location-Upazila he can create new Upazila and edit Upazila.



The screenshot shows a web form titled "Create a New Upazila" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Upazila Name English", "Name Bangla" with a placeholder "Enter Upazila Name Bangla", "GEO Code" with a placeholder "Enter District GEO Code", "District" with a placeholder "Select District", and "Website URL" with a placeholder "Enter Website URL". There is a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button.

29. In Location-Union he can create new union and edit union.



The screenshot shows a web form titled "Create a New Union" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Union Name English", "Name Bangla" with a placeholder "Enter Union Name Bangla", "Upazila" with a placeholder "Select Upazila", "GEO Code" with a placeholder "Enter District GEO Code", and "Website URL" with a placeholder "Enter Website URL". There is a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: a yellow "Close" button and a blue "Create" button, with a red box highlighting the "Create" button.

30. In Location-Ward he can create a new union and edit ward.




The screenshot shows a form titled "Create a New Ward" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter ward Name English", "Name Bangla" with a placeholder "Enter Ward Name Bangla", and a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

31. In Location-Village, can create a new village and edit village.



The screenshot shows a form titled "Create a New Village" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Village Name English", "Name Bangla" with a placeholder "Enter Village Name Bangla", "Union" with a placeholder "Select Union", and a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

32. In Location-Mouza he can create a new mouza, edit, remove mouza.



The screenshot shows a form titled "Create a New Mouza" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Mouza Name English", "Name Bangla" with a placeholder "Enter Mouza Name Bangla", "Division" with a dropdown menu "Select Division", "District" with a dropdown menu "Select District", "Upazila" with a dropdown menu "Select Upazila", "Union" with a dropdown menu "Select Union", and a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.

33. In Location-Pourashava can create a new POURASHAVA, edit, remove POURASHAVA.



The screenshot shows a form titled "Create a New Pourashava" with a close button (X) in the top right corner. The form contains the following fields: "Name English" with a placeholder "Enter Pourashava Name English", "Name Bangla" with a placeholder "Enter Pourashava Name Bangla", "Division" with a dropdown menu "Select Division", "District" with a dropdown menu "Select District", "Upazila" with a dropdown menu "Select Upazila", and a "Status" toggle switch that is currently turned on. At the bottom right, there are two buttons: "Close" and "Create". The "Create" button is highlighted with a red rectangular box, and a red arrow points to it from the left.


34. In Location-City Corporation he can create a new city corporation and edit City Corporation.

35. In Report Section from -the Sample module can get a Sample Wise Report by Filtering with Sample ID, Water Point ID, Caretaker Name, Caretaker Phone Number, and Caretaker NID.




36. System admin can also view and download sample details.


37. In the Lot Wise Report can get the report view Lot wise. System admin He can also download those reports.




38. System admin can get the qualitative report by filtering with different criteria. Such as Date, District, City Corporation, Upazila, Pourashava, Union, Ward, Village, Source of Fund, Project, Type of Water Point, Tender ID, Sample Status, Payment Status, Report Headline.




39. System admin can get the quantitative report by filtering with different criteria. Such as Date, District, City Corporation, Upazila, Pourashava, Union, Ward, Village, Source of Fund, Project, Type of Water Point, Tender ID, Sample Status, Payment Status, Condition Such as Parameter Name, Parameter Value, Parameter Sign. Here the Parameter part is mandatory to see the reports.




40. From the Report Verification module, can enter a Report Memo number and verify. System admin can also download the report.



System admin can download and print the pdf. In the report he will get a QR Code at the bottom of the report. Admin can scan the QR code through mobile and get the report directly in the mobile device.



41. In Management Authority - Admin can create and edit Management Admin.



42. To assign an employee for Management Admin need to enter necessary information

Create Management Admin

Name *
Enter Management Admin Name

Phone No. *
Enter Phone Number

Email
Enter Email

Designation *
Select Designation

Office Of The *
Enter Name of The Office

Status *

Address *

Close Create

43. Select the employee and the start date. At last have to click assign. System admin will get a success message.

WQTRMS

System Admin
Tenggi Zonal Laboratory

Assign New Employee for Management User

User Name
M Admin 2

Employee *
Select Employee

Start Date *
2020/07/2020

Assign Cancel

44. To create a management user , enter the necessary details and click to create.

WQTRMS

System Admin
Tenggi Zonal Laboratory

Management User

View: 50 entries

S#	Chair Name	Information	Designation	Office Of The	Address	Role	Employee	Action
1	M Admin 1	Administrasi	Lab Analisa			Management User	Lab Analisa	Assign
2	M User 2	Administrasi	Substansi	M Admin	270 Elster Garia Road	Management User	Substansi	Assign

Showing 1 to 2 of 2 entries

Create Assign

45. System admin can also see- active or inactive users, remove user, edit management user, and assign an employee for that management user role.

To create a management user System admin has to enter the necessary details and click to create.

Create Management User

Name *
Enter Name

Phone No. *
Enter Phone Number

Email *
Enter Email

Designation *
Select Designation

Office Of The *
Enter Name of The Office

Status *

Address *

Close Create

46. System admin has to assign an employee for Management User. For that he has to click assign.

The screenshot shows a form titled "Create Management Admin" with a close button in the top right corner. The form contains the following fields and controls:

- Name ***: A text input field with the placeholder "Enter Management Admin Name".
- Phone No. ***: A text input field with the placeholder "Enter Phone Number".
- Email**: A text input field with the placeholder "Enter Email".
- Designation ***: A dropdown menu with the placeholder "Select Designation".
- Office Of The ***: A text input field with the placeholder "Enter Name of The Office".
- Status ***: A toggle switch currently turned on.
- Address ***: A large text area for entering the address.
- At the bottom right, there are two buttons: "Close" (orange) and "Create" (blue).

The screenshot shows the "Assign New Employee for Management User" page within the WQTRMS system. The page has a sidebar menu on the left with the following items: Dashboard, Role Manage, Employee Manage, Settings, Lab, Location, Report, Report Verification, and Management Authority. The main content area contains the following fields and controls:

- User Name**: A text input field with the value "M.Admin.2".
- Employee ***: A dropdown menu with the placeholder "Select Employee".
- Start Date ***: A date input field with the value "05/07/2024".
- At the bottom, there are two buttons: "Assign" (red) and "Cancel" (orange).

1.2 Operation Manual of Management Admin


◆ To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login>

The screenshot shows the login page for the "Water Quality Testing & Report Management System". The page features the system logo and the text "Technical Partner: CREATIVE". On the right side, there is a "Sign In" section with the following elements:


- EMAIL**: A text input field with the value "m_admin@gmail.com".
- PASSWORD**: A text input field with masked characters "*****".
- Sign In**: A red button.
- Remember Me**: A checkbox.
- Forgot Password**: A link.
- Download Now**: A link with a QR code and the text "Available on the App Store".

The screenshot shows the "Send to Chemist" dashboard within the WQTRMS system. The sidebar menu on the left includes: Dashboard, Report, Report Verification, and Management Authority. The main content area displays a large number "1" with a red arrow pointing to it, and a "Send" button below it. The user's name "Md.Meris Man" is visible in the top right corner.


From the dashboard management admin will be able to see the management users.




From the Report Module he can see the Lot Wise Report.



From the Report Verification module, he can enter a Report Memo number and Verify. He can also download the report.



From Management Authority, Management admin can create management users



Sixth session

Discussion on Transaction Data Entry System
Customer Registration Entry, Application Management
Sample Collection Form, Bill Management,
Lab Management, Notification Management
Payment Gateway Integration
Development of Mobile Base Application for Data Entry
Responsibilities of System SAE and EE

Sixth session

6.1 Topic:

- ◆ Discussion on Transaction Data Entry System
- ◆ Customer Registration Entry
- ◆ Application Management
- ◆ Sample Collection Form
- ◆ Bill Management
- ◆ Lab Management
- ◆ Notification Management
- ◆ Payment Gateway Integration
- ◆ Development of Mobile Base Application for Data Entry
- ◆ Responsibilities of System SAE and EE

6.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can explain Transaction Data Entry System (Registration, Application Management, Sample collection, Bill Management, LAB Management, Notification, Payment Gateway, Mobile Application)
- ◆ Can understand the responsibilities of System SAE and EE

6.3 Duration: 1.30hr

6.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

6.5 Procedures for Conducting Training

SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none"> ● The name of the subject to be discussed must be presented in writing. ● What the participants want to learn from this discussion should be known through questions. ● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none"> ● The topic should be presented through questions. ● Participants should be encouraged to participate in discussions. ● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none"> ● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none"> ● To assess how much the participants know the subject matter by presenting and answering questions.

6.6.1 Transactional Data Entry

Transactional data relates to the transactions of the organization and includes data that is captured, for example, when a product is sold or purchased. Master data is referred to in different transactions, and examples are customer, product, or supplier data.

1. Customer Registration Entry
2. Application Management
3. Sample Collection
4. Responsibilities of Sub Assistant Engineer (SAE) and Executive (EE)
5. Bill Management
6. Lab Management
7. Notification Management
8. Payment Gateway Integration
9. Development of Mobile Base Application for Data Entry



Training Manual on Water Quality Testing & Report Management (WQTRM) Software



MANUAL DEVELOPED BY:
ENGR. MD. JAHIDUL HASAN, CONSULTANT

FEBRUARY, 2024

Office of the Executive Engineer
Training Department
HRD Building
36-37, Bir Uttam AK Khandkar Road, Mohakhali, Dhaka-1212

ম্যানুয়াল সংকলন

প্রকৌশলী মোঃ জাহিদুল হাসান, কনসাল্ট্যান্ট

মুদ্রণ সহায়তা:

জরুরী ভিত্তিতে রোহিঙ্গা সংকট মোকাবেলায় মাল্টি সেক্টর প্রকল্প (ইএমসিআরপি)

সমন্বয়ক:

প্রশিক্ষণ বিভাগ, জনস্বাস্থ্য প্রকৌশল অধিদপ্তর, মহাখালী, ঢাকা-১২১২

প্রকাশকাল:

নভেম্বর, ২০২৩

ডিজাইন ও মুদ্রণ:

জরুরী ভিত্তিতে রোহিঙ্গা সংকট মোকাবেলায় মাল্টি সেক্টর প্রকল্প (ইএমসিআরপি)

কপি স্বত্বাধিকারী:

জনস্বাস্থ্য প্রকৌশল অধিদপ্তর



Office of the Executive Engineer
Training Department
HRD Building

36-37, Bir Uttam AK Khandkar Road, Mohakhali, Dhaka-1212

1. Customer Registration Entry: All customers who have created online accounts need to provide customer registration information, which is used to capture customer profiles.

2. Application Management:

Create an application module for receiving online application

Create user and application ID

Deliver notification through SMS and e-mail after successful submission of the application

Notification for payment (if pending) status through application code.

3. Sample Collection Form

3. Operation Manual of Sub Assistant Engineer (SAE)

1. SAE can create a water point by clicking on the <https://wqtrms.dphe.gov.bd/water-point> or URL: <https://wqtrms.dphe.gov.bd/login>




2. SAE can create 2 types of Water Point. One is new and another is for existing water point by clicking the “Existence Water Point”, by entering required data




3. Year, Serial Number will be auto generated by clicking create button in Water Point.




4. The Water Point ID is generated as per the Year, District and Upazila Code and Serial Number. By clicking the Create button we can create the Water Point.
5. SAE has to Assign Mechanics to Collect the Water Sample. Can assign two types of mechanic. One is upazila mechanic and another is Lab.
6. In Upazila Mechanic, SAE can choose a mechanic from upazila by dropdown. By clicking assign he can assign that mechanic to collect the sample
7. By selecting Lab, it will be forwarded to lab for assigning sample collector. Similarly have to click Assign Button.




9. In Sample Progress, SAE can see all the Samples. He can view, edit or delete the sample.



10. SAE will get a Sample Collection Form.



11. Information of SAMPLE Area, SAE has to Select District, Upazila, Union, Ward, Village, Mouza, Poura shavas, and City Corporations by Drop-down. Additionally, he has to write Other's Address Name, Other Area Information. By clicking Update Button the Information of SAMPLE Area will be updated.



12. From Caretaker Information of the Water Point, SAE has to enter the Caretaker Name, Caretaker Mobile Number, and Caretaker NID. By clicking Update Button the Caretaker Information of the Water Point will be updated.

[2] Caretaker Information of the Water Point

Caretaker Name *
Enter Household Chief/Caretakers Name

Caretaker Mobile Number *
01 XXX XXX XXX

Caretaker NID
Enter NID Number

Update

X Now Updated

13. From Information About Water Point, SAE will choose Type of Water Point, Ownership of Water Point, Area and Land Type, Type of Water Source, Source of Fund, Project Name, Date of Installation by Drop-down.

14. Additionally, one has to enter manually Installation Depth, Water Source Depth, Filter Location, Filter Length, Distance of Latrine, and Water Point Treatment Unit. There are also two dropdown fields Type of Latrine, Platform of Condition.

15. By clicking Update Button the Information about Water Point will be updated.

16. Sample Information will be updated by the mechanics.

[4] Sample Information

Sample Collection Date
dd--yy--yyyy

Sample Collected By
Mechanic Gazipur 1

Bottle Type
-- Select --

Sample Category
-- Select --

Sample Source
-- Select --

Carried by Sample Box
No

Preservative Used
-- Select --

Latitude


Longitude

Update


X Now Updated

17. From Contractor's Information, SAE has to Choose Tender List by drop-down. According to the Tender the Name, Mobile Number, Package will be auto generated.


18. The Picture Field will be filled by Mechanic during the Sample Collection.




19. In Sample Parameter Assign Field SAE has to choose which parameter they need to send for the test. He can choose the parameter by clicking the check box.




20. After updating all these fields SAE will be able to submit the Sample Application Form. By clicking Final Submit his Sample Application will be submitted.



21. SAE will get a Success message and In Sample Progress he can view, edit, and forward and delete the Water Point/Sample.



22. For further process, SAE has to forward it to the Executive Engineer Panel. By Clicking the Forward Button, he will be able to send it to the EE panel.




23 After Forwarding SAE will get a success message that the sample forward to the Executive Engineer panel.


4. Operation Manual of EXECUTIVE ENGINEER (EE)



1. To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login>




2. The Executive Engineer will get a notification; that the sample is forwarded to the
3. EE will get the Forwarded Sample in Sample Progress. He can see the Sample Details and also can backward it to the SAE panel.






4. To generate lot, EE has to select Tender ID, Upazila, Sample ID, and Sample Collect by, Project from drop-down. If the Tender ID is null, he can choose the Null Tender Checkbox.

5. There will be also a list of samples that was forwarded to SAE. EE can select multiple water point or any particular water point to generate lot by clicking checkbox.



6. After generating a lot, he will get a Success Message.




10. Other details will be auto selected. To complete the payment, he has to click the Submit Payment button.

11. EE will get a message after successful payment.



12. At lasted has to forward the lot to the lab receiver for further processing. He can Select Lab by dropdown and forward it.

13. Executive Engineer will get a success message after forwarding it.



Seventh session

Discussion on Transaction Data Entry System
Customer Registration Entry, Application Management
Sample Collection Form, Bill Management,
Lab Management, Notification Management
Payment Gateway Integration
Development of Mobile Base Application for Data Entry
Responsibilities of System SAE and EE

Seventh session

7.1 Topic:

- ◆ Lab Management System (LIMS)
- ◆ Responsibilities of LAB Technician, Chemist, LAB In-Charge
- ◆ Notification Management System
- ◆ Mobile Application

7.2 Objectives of the lesson: Participants at the end of this session -

- ◆ Can learn about Lab Management Systems (LIMS)
- ◆ Can Learn about the Responsibilities of LAB Technician, Chemist, LAB In-Charge
- ◆ Can learn about Notification Management System
- ◆ Can learn about Mobile Application

7.3 Duration: 1.00hr

7.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

7.5 Procedures for Conducting Training


SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none">● The name of the subject to be discussed must be presented in writing.● What the participants want to learn from this discussion should be known through questions.● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none">● The topic should be presented through questions.● Participants should be encouraged to participate in discussions.● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none">● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none">● To assess how much the participants know the subject matter by presenting and answering questions.

7.6.1 Lab Management System (LIMS):

A Laboratory Information Management System (LIMS) is software that allows you to effectively manage samples and associated data. By using a LIMS, your lab can automate workflows, integrate instruments, and manage samples and associated information.

The main purpose of a LIMS is to improve lab efficiency, quality control, and accuracy by reducing manual operations and generating LAB reports.


- ◆ Lab Profile/Lab Office Setup
- ◆ Lab user profile
- ◆ Sample Collection Scheduler
- ◆ Sample and Test type
- ◆ Sample test parameter builder
- ◆ Test result builder
- ◆ Attach inspected data or sample image
- ◆ Test report builder




LIS Flow Chart

Operation Manual of LAB Admin


1. To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login>




2. In Lab-List Admin can Edit Lab Information, View Lab Parameters, View Lab Users, View Lab Details, View Lab Administer Committee, Remove Lab,




3. Lab Admin can edit lab information. Lab Admin can also create a new lab.



4. Lab Admin can view and add parameters.



5. Lab Admin can view lab users, edit lab users, and assign employees.



6. Lab Admin can also create a new lab.



Operation Manual of LAB In Charge


To Create Water Point:

1. To Create Water Point, use URL: <https://wqtrms.dphe.gov.bd/login>



2. To log in, please use the provided User Name and Password.

3. Lab In charge can create a water point by "Water Point" module. Alternatively, Lab In charge can directly access the Water Point section by going to the following URL: <https://wqtrms.dphe.gov.bd/water-point>




4. Lab In charge can create three types of Water Point-1. newly Create Water Point 2. for Existence Water Point 3. Manually Water Point.


5. For Create Water Point Lab In charge has to Choose Year, District, Upazila and Serial Number will be auto generated. By clicking, create button he will be able generate Water Point.

6. Lab In charge will get a Success message.


7. For Existence Water Point he has to Choose Water Point Number and click create to create the water point.




8. Water Point, he has entered a 23-digit ID to generate and click the create button. Lab In charge will get a Success message.



Lab In charge will get a Success message.




9. To Assign Sample Collector, he has to click to Assign Sample Collector. Then Lab In charge has to choose to Lab Sample Collector and Assign Sample collector.




10. Lab In charge will get a notification after Collecting the Sample by Lab Sample Collector.




- It will be stored in Sample Progress. Here Lab In charge can view Sample Details, Edit Sample Form and Delete the Sample.
- To edit Lab In charge has to click the second option.




- From Information of Sample Area, Lab In charge has to Select District, Upazila, Union, Ward, Village, Mouza, Pourashavas, and City Corporations by Drop-down. Additionally, he has to write Other's Address Name, Other Area Information. By clicking the Update Button, the Information of Sample Area will be updated.



- In Caretaker Information of the Water Point, Lab In charge has to enter the Caretaker Name, Caretaker Mobile Number, and Caretaker NID. By clicking Update Button the Caretaker Information of the Water Point will be updated.



- In Information about Water Point, Lab In charge will choose Type of Water Point, Ownership of Water Point, Area and Land Type, Type of Water Source, Source of Fund, Project Name, Date of Installation by Drop-down.
- Additionally, one has to enter Installation Depth, Water Source Depth, Filter Location, Filter Length, Distance of Latrine, and Water Point Treatment Unit. There are also two dropdown fields Type of Latrine, Platform of Condition. By clicking the Update Button, the Information about Water Point will be updated.



17. Sample Information will be updated by the Lab Sample Collector.

[4] Sample Information

Sample Collection Date: dd-yyy-yyyy

Sample Collected By: Mechanic Gasipur 1

Bottle Type: -- Select --

Sample Category: -- Select --

Sample Source: -- Select --

Carried by Sample Box: No

Preservative Used: -- Select --

Latitude:

Longitude:

x Now Updated

[5] Contractor's Information

Tender List: -- Select --


Name:

Mobile Number:

Package No:

x Now Updated

[6] Picture




[7] Sample Parameter Assign

Select Parameters From Below Panel




<input type="checkbox"/> Aluminum(Al)	<input type="checkbox"/> Ammonia
<input type="checkbox"/> Arsenic(As)	<input type="checkbox"/> Barium(Ba)
<input type="checkbox"/> Boron	<input type="checkbox"/> Calcium(Ca)
<input type="checkbox"/> Chemical Oxygen Demand(COD)	<input type="checkbox"/> Chloride
<input type="checkbox"/> Chlorinated(2, 4, 6 Trichlorophenol)	<input type="checkbox"/> Coliform(Faecal)
<input type="checkbox"/> Coliform(Total)	<input type="checkbox"/> Colour
<input type="checkbox"/> Iron(Fe)	<input type="checkbox"/> Manganese(Mn)
<input type="checkbox"/> PH	<input type="checkbox"/> Total Dissolved Solid(TDS)
<input type="checkbox"/> Total Hardness(as CaCO3)	<input type="checkbox"/> Turbidity

x Now Updated



21. After updating all these fields Lab In charge will be able to submit the Sample Application Form. By clicking Final Submit his Sample Application will be submitted.




Approve Process:





4. To confirm the forward, he has to click the forward button.



Operation Manual of LAB Receiver




2. The Lab Receiver will get a notification; that the lot is forwarded to them.
3. He can See the Payment Status, Received Status, and Backward to Executive Engineer, View the Lot Samples.
4. To receive the lot Lab Receiver has to click the second option.
5. To receive the lot Lab Receiver has to select the Status by drop-down and make a Comment and finally submit by clicking Submit Receive.




Operation Manual of LAB Analyzer


1. To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login>



2. The Lab Analyzer will get a notification that the lot is forwarded to them.




3. Lab Analyzer can see the Payment Status, Received Status and the Lot Samples.




4. In Sample Progress, he will get the Samples. He can view Sample Details and Sample Test.


5. To test the sample, he has to take the View Sample Test.




6. Analyzer will get the Parameters that are needed to be tested. In Each Parameter, he/ she has to select the Start Date, End Date, Time, Analysis Procedure, Analysis Method, and Tested By in dropdown. It has to manually enter the Concentration, LOQ and Comment.



7. At last, he/ she has to Save the Result by Result Save option and He/she will get a Success message after submitting the result.



- 8. After testing the parameters, he will get a Comment option and a Test Result Submit Button.
- 9. Lab Analyzer will get a success message after submitting the test result.




Operation Manual of Chemist (Junior Chemist, Chemist, Senior Chemist, Chief Chemist)


1. To access WQTRMS, Use URL: <https://wqtrms.dphe.gov.bd/login>




2. The Chemist Panel will get a notification that the lot is forwarded to them.



3. Chemist will get the Sample in Sample Progress. Here he can View Sample Details, Check Sample Test Status and Click for Approval.



4. To approve the sample, he has to click for approval.



5. Here Chemist will get a Status and Date to select and a Comment box to comment.



6. Chemist has to click the Update button to update the Status.


7. Chemist will get a Success Message.



Note: Approval process is same for all chemist. It will sequentially go from the Junior Chemist, Chemist, Senior Chemist, Chief Chemist. Only one Chemist can publish the report according to the Lab Policy.

Report Publish:

- To publish the report Chemist will go to the Sample Forward module where he can click to publish the report.



- Chemist will be able to see the Approval from the Lab Authority, he has to select Status and Date. In Comment section he has to write the Recommendation. By clicking Update Button he will be able to generate the report.



- Chemist will get a success message after publishing.



- Chemist will get the published report in the Sample-Publish module.



Eight session

**Payment Gateway
Mobile Base Application
Notification Management System
Dynamic Report
Dynamic Dashboard**

Eight session

8.1 Topic:

- ◆ Payment Gateway
- ◆ Mobile Base Application
- ◆ Notification Management System
- ◆ Dynamic Report
- ◆ Dynamic Dashboard

8.2 Objectives of the lesson: Participants can learn the followings at the end of this session –

- ◆ Payment Gateway
- ◆ Mobile Base Application
- ◆ Notification Management System
- ◆ Can learn about Dynamic Report
- ◆ Can learn about Dynamic Dashboard
- ◆ Can learn about the Challenges of the new system

8.3 Duration: 1.00hr

8.4 Training Materials: Whiteboard, whiteboard marker, multimedia screen, duster, and manual.

8.5 Procedures for Conducting Training

SUBJECT	ACTIVITIES
Name and purpose of discussion	<ul style="list-style-type: none"> ● The name of the subject to be discussed must be presented in writing. ● What the participants want to learn from this discussion should be known through questions. ● Pre-written objectives should be presented and participants' expectations should be compared with them.
Discussion	<ul style="list-style-type: none"> ● The topic should be presented through questions. ● Participants should be encouraged to participate in discussions. ● Care must be taken whether questions are presented systematically or not.
Summary of topics discussed	<ul style="list-style-type: none"> ● At the end of the discussion, the entire subject should be discussed.
Assessment	<ul style="list-style-type: none"> ● To assess how much the participants know the subject matter by presenting and answering questions.

8.6.1 Dynamic Report System:

- ◆ The proposed system has a dynamic report engine that will generate reports based on
- ◆ dynamic data as per clients' requirement
- ◆ The system has drag and drop report builder engine to facilitate any kind of user/system report
- ◆ The proposed system has an option to generate various reports for analysis
- ◆ The proposed system has the functionality to correct or modify the generated report
- ◆ Basic reports for different types/levels of user containing no. of sample collected, no. of sample tested, office/lab wise sample collection, office/lab wise test result, no. of emergency report, etc.
- ◆ The proposed system has the provision for generating different analytical reports based on
- ◆ assessments with graphical representations like charts, pies, lines, etc.
- ◆ List of services with a price/ center with address/Users with role
- ◆ Division/District/Upazila-wise client list
- ◆ Date to Date bill collection reports Detail/SUM
- ◆ Date to Date Division/District/Upazila wise test report detail/SUM

8.6.2 Dynamic Dashboard:

The dashboard will contain statistical data on different features of the system in different interactive formats like indicator-based graphical charts, tables with figure/summary data, etc., which are easy to understand at a glance.

The Dashboard can present dynamic data of sample testing such as the number of enrolled, number of completed, number of not collected, results distribution among the applicant's types, etc. It can present a zoning map of metro areas, municipalities, growth centers, rural areas, and the whole country based on clustering and Map for the risk level pattern of the water points. The system will display an office-based dashboard like

- Dashboard for Operating Office-Upazila/Union, District Office, Division Department/HQ

8.6.3 Challenges of the proposed system:



8.6.4 Notification Management System:

Notification Management helps you ensure the accurate delivery of important information to the right person at the right time, without depending on employees to take action. By using Notification Management, you can notify customers immediately regarding order confirmations, delayed deliveries, or changing specifications.

- The system has provision for Notification through SMS, e-mail and system
- The system or admin user will be able to configure notification media/content
- Able to configure notification template
- The proposed system has a provision to keep a log of every notification generated or sent through notification media
- Notification for citizen application status, history progress, etc

8.6.5 Payment Gateway Integration:

The gateway is the beginning and end of the transaction. The customer will enter their credit card information and receive an approval or denial of the transaction. The payment processor moves the information between the customer's bank and the merchant bank. Every transaction processed online needs both.

To integrate payment gateway for online payment collection and also payment settlement facilities for respective offices with notification and detailed reports. Notify after payment through SMS and e-mail.



8.6.6 Mobile Application:






A mobile app (or mobile application) is a software application developed specifically for use on small, wireless computing devices, such as smartphones and tablets, rather than desktop or laptop computers.




Uses of Mobile Application in WQTRM System

- Development of Mobile application for sample collection, lab result entry, view and download report
- The mobile app can download the citizen application to mobile local storage while connected to the server with the internet.
- The mobile app can accept input of sample test result data offline and synchronize the data with the server when the mobile is connected to the internet later.
- The mobile app will accept sample test result data either online or offline only if GPS is on. The mobile app can process water point GPS information from the picture taken during ample testing.

10. Mobile Apps for Mechanics & Lab Sample Collector


1. To access WQTRMS Mobile Apps, please open our provided APK or download the Mobile Application from the web link <https://wqtrms.dphe.gov.bd/login>





Operation Manual of Report Management System


1. To access WQTRMS, use URL: <https://wqtrms.dphe.gov.bd/login>



2. In the Dashboard authority can see the Sample Report, Lot Report, Qualitative and Quantitative Report.



3. In the Report Section from the Sample module, he can get Sample Wise Report by Filtering with Sample ID, Water Point ID, Caretaker Name, Caretaker Phone Number, and Caretaker NID.




4. In the Lot Wise Report, he will get the report Lot wise and can download the report.



5. Authority can download and print the report.



6. Authority can get the qualitative report by filtering with different criteria. Such as Date, District City Corporation, Upazila, Pourashava, Union, Ward, Village, Source of Fund, Project, Type of Water Point, Tender ID, Sample Status, Payment Status, Report Headline.



7. Authority can Export the qualitative report and it will be downloaded in excel file



8. Authority can select item for PDF and select PDF layout.



9. Then authority will be able to print and download the report.



10. Authority can download report in PDF with QR code and also in excel file as per parameter.



11. From the Report Verification module, he can enter a Report Memo number and verify. He can also download the report.



12. Authority can also download the report and print the pdf with QR Code.



13. Authority can also download the report and print the pdf with QR Code.

