

## BRiCM LABS® EZ Dengue Ag

### FEATURE

- Quick result (5-10 min)
- Easy to use
- Can be done at home
- Long self-life
- No special requirement for storage or transportation
- No specialized equipment is needed

### SPECIFICATIONS

- Reliable detection of all four serotypes (DEN1, 2, 3, 4)
- Sensitivity and specificity: 100% against gold standard
- Limit of detection: 0.25 ng/mL
- Contains distinct C (control) line to validate the result
- No cross-reactivity with Chikungunya, Covid-19 and ZIKV
- Self-life: 24 months

### INTENDED USE

- For the detection of Dengue virus from human whole blood, serum, and plasma

### PACK SIZE

Each box contains

- 25 pouches with rapid cassette
- 25 graduated plastic dropper
- 25 sterile spring loaded lancet
- 3ml sample dilution
- Manual

### STORAGE

- Store at 4°C - 50°C

### CAUTIONS

- Do not use the same strip twice
- Use immediately after opening the pouch
- Do not use it if the pouch is open
- Only use as recommended by the instruction

